


biomarkers for early-stage alzheimer's disease

By: Dr Jacqueline Gosink, EUROIMMUN
AG, Luebeck, Germany



the middle stage, patients struggle to manage everyday life; they become confused and suffer from hallucinations amongst other things. By the final stage, patients are increasingly frail. They need constant support and care and have difficulties in eating and sometimes in swallowing. On average, the life expectancy after the onset of symptoms is seven to ten years.

PATHOLOGICAL CHARACTERISTICS

Alzheimer's disease is characterised pathologically by the formation of deposits in the neuronal cell body and outside the nerve cell ends (see *figure 1*). The intracellular deposits (neurofibrillary tangles) consist of hyperphosphorylated tau proteins (P-tau), which form tangled fibre strands. The extracellular deposits (neuritic plaques) contain predominantly the peptides A β 1-40 and A β 1-42, which are breakdown products of the membrane-bound amyloid precursor protein. The physiological function of A β 1-40 and A β 1-42 peptides is not yet fully understood, but they are presumed to play a major role in signal transmission in nerve cells.

DIAGNOSIS

Definitive diagnosis of Alzheimer's disease can only be established post mortem. In autopsy the neuropathological changes, namely plaques and neurofibrillary tangles, are visible in the brain of the deceased patient. Tentative *in vivo* diagnosis of Alzheimer's disease ('probable Alzheimer's disease') is based primarily on the clinical sign of memory loss and the exclusion of possible reversible causes. Imaging techniques such as MRT, SPECT or PET (amyloid detection) can be used to support early and differential diagnostics, for example exclusion of endocrinopathies and electrolyte disorders. The results from imaging dementia diagnostics are assessed together with other available diagnostic information, including analyses of CSF.

CSF BIOMARKERS

The diagnosis of Alzheimer's disease in the early and pre-symptomatic stages requires reliable, quantifiable CSF biomarkers, for example soluble A β 1-42 and tau proteins. The concentrations of these analytes in the CSF reflect the neuropathological changes in the brain. There are currently no blood markers available that show the same clinical value as the CSF markers. Patients with Alzheimer's disease show a significantly decreased level of A β 1-42 that is already detectable 5-10 years before the start of cognitive changes. The concentrations of total tau and P-tau, on the other hand, increase when patients show advanced neurodegeneration and cognitive impairment. It is thus possible to discriminate Alzheimer's patients from healthy persons by means of CSF markers.

In contrast to A β 1-42, the level of A β 1-40 remains unchanged in Alzheimer's patients. Measuring A β 1-40 together with A β 1-42 enables calculation of the A β 1-42 to 1-40 ratio, which helps to increase the efficiency of early diagnostics (see *figures 2 and 3*). A ratio of under 0.1 indicates amyloid pathology. This ratio might further help to discriminate Alzheimer's disease from other diseases, such as vascular dementia. →


Alzheimer's disease is a major global challenge of our time, with its prevalence predicted to rise dramatically as population's age. No cure is available, and early diagnosis is crucial for managing the disease. Since clinical diagnosis is difficult, especially in the initial stages, analysis of biomarkers in patients' cerebrospinal fluid (CSF) is increasingly used as a diagnostic support tool. Analysis of the beta-amyloid (A β) peptides 1-42 and 1-40 together with total tau and phosphorylated tau (P-tau) can significantly aid early diagnosis of Alzheimer's disease. These neurochemical markers can be measured precisely, reproducibly and independently of matrix effects using a new generation of ELISAs based on well-characterised capture antibodies.

ALZHEIMER'S DISEASE


Alzheimer's disease was first described in 1906, and is with 60-70% the most common cause of dementia in old age. The risk for developing Alzheimer's disease doubles for around every five years after age 65, with 30% of persons over 90 suffering from the disease. In contrast to the age-dependent, sporadic form of Alzheimer's disease, the familial, genetically caused form can occur in young adults from 30 years of age.

The disease course is divided into three consecutive phases; the pre-clinical stage, the mild cognitive impairment (MCI) stage and the dementia stage. Patients usually attend a neurologist initially because of learning difficulties and short-term memory impairment. Cognitive function then declines progressively. In


▼ **FIGURE 1:** Neuron showing pathological characteristics of Alzheimer's disease




▼ **FIGURE 2:** Determination of Aβ concentrations and ratio calculation




▼ **FIGURE 3:** Differential diagnostics using Aβ ratio



▼ **FIGURE 4:** Epitopes of Aβ antibodies



▼ **FIGURE 5:** Principle of the ELISA technology



and automatable, allowing CSF diagnostics to be easily integrated into the automated routine operations of a diagnostic laboratory.

CLINICAL EVALUATION

Samples from clinically characterised patients with Alzheimer's disease (85) or vascular dementia (10) and control subjects (18) were analysed using the Aβ 1-42 and 1-40 ELISAs (see figure 3). 82 of the Alzheimer's patients yielded an Aβ 1-42 value of below 400 pg/ml and a Aβ ratio of below 0.1, in line with the clinical diagnosis. In contrast, patients with vascular dementia exhibited Aβ 1-42 values that were on average twice as high (422 pg/ml) as those of Alzheimer's patients (210 pg/ml) and an Aβ ratio of 0.10 to 0.17. The control group showed a mean Aβ 1-42 value of 543 pg/ml and a ratio of between 0.12 and 0.17. These data demonstrate the usefulness of determining Aβ 1-42, Aβ 1-40 and the ratio of the two analytes in the diagnosis and differentiation of Alzheimer's disease.

PERSPECTIVES

With the current worldwide prevalence of 36 million Alzheimer's cases expected to double by 2030, the need to reliably identify patients and individuals at risk is critical. The advent of CSF assays for biomarkers, such as beta amyloid peptides and tau proteins, represents a significant step forward in tackling the disease. Molecular genetic risk determination, for example the detection of Alzheimer's-associated apolipoprotein E alleles by DNA microarray, may also prove to be helpful in early diagnostics. Neurochemical and genetic indicators are likely to remain at the forefront of Alzheimer's research as scientists hunt for strategies that allow ever-earlier diagnosis. Predicting Alzheimer's even before it manifests would offer the best chance of preventing its devastating consequences. 